LANAFORM

HEALTH

ABPM-100

ARM BLOOD
PRESSURE
MONITOR

4+30
MEMORIES

WHO RISK
INDICATOR

IRREGULAR HEART
RHYTHM ALERT

MEDICAL
DEVICE

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NADLAKTNI MERILNIK
KRVNEGA TLAKA

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EN Instruction manual

📷 *The photographs and other images of the product in this manual and on the packaging are as accurate as possible, but do not claim to be a perfect reproduction of the product.*

Thank you for choosing the LANAFORM ABPM-100 arm blood pressure monitor. The arm blood pressure monitor is designed to be a non-invasive measuring process and for monitoring adults' blood pressure. It allows you to measure your blood pressure quickly and easily, register the values and display changes in values. It warns you if you have an irregular heart beat.

TITLE AND DESCRIPTION OF DISPLAY SYMBOLS ON THE INSERT DESCRIBING THE PRODUCT'S FEATURES

Definition of the symbol	Identification	Information provided by the manufacturer
Manufacturer's details		Manufacturer
Batch number		Production batch number
BF device		The cuff is type BF applied part
Instructions for use		Follow the instructions for use The sign background colour: blue. The sign graphical symbol: white
Please note: see instructions for use		Please note: consult the enclosed documents. Please read these instructions before using the device. Please also follow these instructions carefully when using the device.
Notified body		
European representative		European representative
Electrical and electronic waste (WEEE)		The device should not be disposed of with other household waste at the end of its serviceable life.
Date of manufacture		
Serial number		

Keep dry

IP20

IP20

The first characteristic numerical symbol for "Degrees of protection against access to hazardous parts and against solid foreign objects". The second characteristic numerical symbol for "Degrees of protection against ingress of water"

▲ PLEASE READ ALL THE INSTRUCTIONS BEFORE USING YOUR BLOOD-PRESSURE MONITOR, PARTICULARLY THESE FEW BASIC SAFETY INSTRUCTIONS

- When using this product, please make sure you comply with all the items listed below. Any deviation from these instructions may cause injury or affect the accuracy of the readings.
- Always measure your blood pressure at the same time of day so that the values are comparable.
- Rest for about 5 minutes before taking any measurements.
- When you have to take several measurements on an individual, wait for 5 minutes between each measurement.
- Do not eat, drink, smoke or take part in physical activities for at least 30 minutes beforehand.
- Refrain from talking when the monitor is recording measurements in order to guarantee their accuracy.
- Take another measurement if you are at all in doubt about the values you have been measuring.
- The measurements you have taken are purely for information purposes – they are not a substitute for a proper medical examination!
- Notify your doctor of your results since you should never under any circumstances take any medical decisions based on these measurements alone (for example, choosing medicines and their dosages).
- Do not use the blood-pressure monitor on new-borns, pregnant women and patients with pre-eclampsia.
- It is inappropriate for people with serious arrhythmia to use this Electronic Sphygmomanometer.

- Do not use the cuff other than supplied by the manufacturer, otherwise it may bring biocompatible hazard and might result in measurement error.
- If the blood circulation is restricted in an arm due to chronic or severe vascular diseases (including vasoconstriction), the accuracy of the measuring process on the wrist is limited. In this case, use an upper arm blood pressure monitor.
- Cardio-vascular disease can cause incorrect or inaccurate measurements. This is also true for very low blood pressure, diabetes, circulation problems, irregular heartbeats and chills, fever or shaking.
- The monitor must not be used at the same time as any high-frequency surgical device.
- Only use the device on people whose arm perimeter is the same as that stated for the device.
- Please note that the function of the limb concerned may be hampered when inflating.
- Blood circulation must not be blocked for longer than necessary when taking blood pressure. If the device does not work, remove the cuff from the arm.
- Do not take measurements too frequently and avoid sustained pressure in the cuff. This causes blood circulation to be reduced and may cause injury.
- Please do not put the cuff on an arm, the arteries or veins of which are undergoing medical treatment, e.g. where an intravascular access device is present which is intended for intra-vascular treatment or if there is an arteriovenous shunt in place.
- Do not use the cuff on people who have undergone a mastectomy.
- Do not put the cuff over wounds as its use may aggravate them.
- Put the cuff only round the arm. Do not put the cuff over other parts of the body.
- The blood-pressure monitor is battery-operated only.
- Automatic switch-off means that the monitor will go into energy-saving mode when left untouched for one minute.
- The device is designed for use as described in these instructions for use. The manufacturer cannot be held responsible for damage caused by inappropriate or improper use.
- Do not use the cuff other than supplied by the

manufacturer, otherwise it may bring biocompatible hazard and might result in measurement error.

- Please do not share the cuff with other infective person to avoid cross-infection.
- Avoid compression or restriction of the connection tubing during measurement, which may cause inflation error, or harmful injury due to continuous cuff pressure.
- This Electronic Sphygmomanometer is designed for adults and should never be used on infants or young children. Consult your physician or other health care professionals before use on older children.

▲ Storage and maintenance instructions

- The monitor is made up of precision electronic components. The accuracy of the values measured and the serviceable life of the device depend on careful handling:
 - Protect the device from shocks, humidity, dirt, marked fluctuations in temperature and direct sunlight.
 - Do not drop the device.
 - Do not use the device in the vicinity of strong electro-magnetic fields and keep it away from cordless or portable telephone facilities, etc.
- It is recommended that the batteries be removed when the device is not being used for long periods.

▲ Notes on batteries

- If the liquid of the battery cell comes into contact with the skin or eyes, rinse the area affected with water and consult a doctor.

▲ *Ingestion hazard! Young children might swallow the batteries and suffocate. So keep the batteries out of the reach of children!*

- Check the plus (+) et minus (-) polarity symbols.
- If a battery has leaked, put on protective gloves and clean the battery compartment with a dry cloth.
- Protect the batteries from excessive heat.


⚠ Explosion hazard! Do not throw batteries onto the fire.

- Do not attempt to re-charge the batteries.
- Remove the batteries from the battery compartment if the device is not to be used for a lengthy period.
- Use only identical or equivalent batteries.
- Always replace all the batteries at the same time.
- Do not use rechargeable batteries.
- Do not dismantle, open or break the batteries.

▲ Instructions for repair and removal


- Batteries must not be disposed of with household waste.
- Please take flat batteries to the collection points provided for this purpose.
- Do not open the device. Failure to follow these instructions may void the warranty.
- You must not repair or adjust the device yourself. If you have done so, safe operation of the device can no longer be guaranteed.
- Only the customer service department or an authorised operator may carry out any repairs.
- Before filing any claim, however, please check the battery status beforehand and replace them, if necessary.
- In the interests of environmental protection, the device should not be disposed of with other household waste at the end of its serviceable life. It should be disposed of via your country's authorised collection points. The device must be disposed of in accordance with the European WEEE (Waste Electrical and Electronic equipment) Directive on used electrical and electronic devices. If you have any questions on this, please contact the local bodies responsible for the disposal and recycling of these products.

DESCRIPTION OF THE DEVICE

- 1 Screen
- 2 Stop/start button
- 3 Memory button
- 4 Cuff connection socket
- 5 Battery compartment
- 6 Scale for ranking measuring results
- 7 Cuff
- 8 Cuff pipe
- 9 Connection to the cuff

DISPLAY SCREEN

- 10 Systolic pressure
- 11 Diastolic pressure
- 12 Pulse rate measurement
- 13 Pulse symbol
- 14 User memory
- 15 Backup location number/average memory value (AL), morning (A), evening (P)
- 16 Memory function symbol
- 17 Irregular heart beat symbol
- 18 Ranking measuring results
- 19 Deflation
- 20 Time et date

21 Battery change symbol

INSERTING THE BATTERIES

- Remove the battery compartment cover (5) on the left-hand side of the device.
- Insert four 1.5V micro batteries (LR03 alkaline not included in the packaging).
- You must insert the batteries according to the polarity illustrated. Do not use rechargeable batteries.
- Close the battery compartment cover carefully.
- If the battery change symbol (21) appears, no measurements can be taken and you need to change all the batteries.

▲ Removing the batteries

- Flat and fully discharged batteries must be disposed of in special containers or at collection points intended for this purpose or returned to an electrical appliance retailer. You have a legal obligation to dispose of batteries properly.

SET-UP DATE AND TIME

- This menu allows you to adjust the following functions one after the other:
Time → Date → Hour format

▲ *You must adjust the date and time. Without these adjustments, it is impossible to register and read your measurements accurately with the date and time.*

Time format

- Press and hold down the STOP/START button (2) and the memory button (3) at the same time for 3 seconds or insert the batteries.
- Using the memory button (3), select the time format you require and confirm by pressing the STOP/START button (2).

Date

- The year flashes on the screen. Using the memory button (3), select the year you require and confirm by pressing the STOP/START button (2). 
- The month flashes on the screen. Using the memory button (3), select the month you require and confirm by pressing the STOP/START button (2). 
- The day flashes on the screen. Using the memory button (3), select the day you require and confirm by pressing the STOP/START button (2). 


▲ *If the time format has been adjusted to 12 pm, the order of the day and month displayed is reversed.*

Time

- The time flashes on the screen. Using the memory button (3), select the time you require and confirm by pressing the STOP/START button (2). 
- The minutes flash on the screen. Using the memory button (3), select the minutes you require and confirm by pressing STOP/START (2). 
- Once all the dates have been adjusted, the device automatically switches off.

MEASURING BLOOD PRESSURE

Positioning the cuff

- Put the cuff on your bare left arm. The circulation in your arm must not be restricted by clothes that are too tight or anything else.
- Position the cuff on your arm so that the lower edge 

is 1 - 2cm above the crook of your elbow and above the artery. The pipe is directed towards the middle of your palm.

- Now close the free part of the cuff around your arm using the Velcro without making it too tight. You must be able to insert two fingers under the cuff.
- Now connect the cuff's pipe into the cuff's power supply.

▲ The device must only be used with the original cuff. The cuff is intended for an arm circumference of 22 to 30cm.

Adjusting to the correct position

- Rest for about 5 minutes before taking any measurements. Otherwise, this may cause errors.
- Measurements can only be taken when sitting down or prone. Make sure in every case that the cuff is at heart level.
- Make yourself comfortable before taking your blood pressure. Make sure that your back and your arms are well supported on the backrest and armrests. Do not cross your legs. Put your feet flat on the floor.
- In order not to distort the result, it is important to remain calm when taking the measurements and not to talk.

Measuring blood pressure

- Position the cuff as described above and get into the position you want to take the measurements.
- To start the monitor, press the STOP/START button (2). All displays will light up briefly.

- To choose the memory location you want, press the memory function button (3). (1, 2, 3) or (4).
- The results of measuring 4 different people or the morning and evening measurements may be recorded separately in four memories with 30 places.
- The monitor automatically starts the measuring process.
- You can stop the measuring process at any time by pressing the STOP/START button (2). As soon as it recognises a pulse, the Pulse symbol is displayed.

- The results of measuring the systolic pressure (10), diastolic pressure (11) and the pulse (12) are displayed. The irregular heartbeat symbol (if there is one) (17) and the symbol for ranking measuring results (18) flash.
- Er_ is displayed when the measuring process has not been carried out correctly. Read the chapter Error message/Rectifying errors in these instructions for use and start taking measurements again.
- Switch off the monitor by pressing the STOP/START button (2). The measuring result is then recorded in the user memory you have chosen.
- If you forget to switch off the device, it will switch off automatically after approximately 1 minute.
- Wait at least 5 minutes before taking another measurement.

EVALUATING THE RESULTS

Irregular heart beat

- This device is able to identify any irregular heart beat during the measuring process and displays it, if necessary, after the measuring process with the corresponding symbol (17).

- This may indicate arrhythmia. Arrhythmia is a disease characterised by an anomaly in the heart rate due to disturbances in the bio-electrical system. The symptoms (heart beats that are late or early, faster or slower pulse) can be caused particularly by heart disease, by age, by natural predisposition, by a diet that is too rich, by stress or even by lack of sleep. Only a medical consultation can determine arrhythmia. If the symbol is displayed after the measuring process, start again. Make sure that you rest for 5 minutes and do not talk or move during the measuring process. If the symbol often appears, please consult a doctor. Self-diagnosis and self-medication based on the results of the measuring process can be dangerous. You must follow your doctor's instructions.

Ranking measuring results


- The results of the measurement process are ranked and evaluated according to the table below. These values must only be used for information purposes because blood pressure differs from person to person and depending on their age, etc. It is important to consult your doctor regularly. Your doctor will provide you with your own values for normal blood pressure and the level at which blood pressure is regarded as dangerous. The ranking is displayed as is the scale of the device for establishing the range within which blood pressure has been measured. If the systolic and diastolic values are within two different ranges (e.g. systolic within the "normal high" range and diastolic within the "normal" range), the scale markings for the chart will always give the highest range on the device, i.e. "normal high" in this example.

Range of blood pressure values	Systolic (in mmHg)	Diastolic (in mmHg)	Colour displayed
Level 3: acute hypertension	≥ 180	≥ 110	Red
Level 2: intermediate hypertension	160 – 179	100 – 109	Orange
Level 1: mild hypertension	140-159	90 - 99	Yellow
Normal high	130 - 139	85 - 89	Green
Normal	120 - 129	80 - 84	Green
Optimum	<120	<80	Green

Source: WHO, 1999 (World Health Organisation)

Recovering and removing measuring values

- The result of every successful measurement is recorded with the date and time. When there are more than 30 measuring data, the oldest measurement is deleted.
- Press the memory button (3) to display the measuring result. The number of results recorded for the selected user is displayed.

- Press the STOP/START button (2) to select another user and confirm by pressing the memory button (3).
- The average value (AL) (15) of all the values measured and recorded for this user is then displayed.

- If there is no result in the memory of the selected user, the LCD screen displays "0" for blood pressure and for the pulse measurement.

- Press the memory button (3). A7 flashes on the screen. The average value of morning measurements of the last 7 days is displayed (morning: 5 am – 9 am).

- Press the memory button (3). P7 flashes on the screen. The average value of evening measurements of the last 7 days is displayed (evening: 6 pm 8 pm).

- If you press the memory button again (3) again, the latest measurement is displayed on the screen. The irregular heartbeat symbol (if there is one) (17) and the symbol for ranking the measuring results (18) flash.

- If you press the memory button (3) again, you can see your individual measurements.
- To switch the device off, press the STOP/START button (2).
- You can leave the menu at any time by pressing the STOP/START button (2).

Deleting the values measured

- To delete the records of a specific user memory, first select a user memory.
- Hold the memory button (3) down for 3 seconds to delete the results.

ERROR MESSAGE/RECTIFYING ERRORS

Problem	Cause	Solution
The LCD screen displays an abnormal result.	The cuff is in the wrong position or has not been tightened sufficiently.	Apply the cuff correctly and try again.
	The body's posture is not correct during the test	Re-read sections "Adjusting to the correct position" in the instructions and try again.
Arm or body movement, nervousness, excitement or the act of talking during the measuring process	Heart rate is irregular (arrhythmia)	Take another measurement and follow the positioning instructions and instructions for use. The blood pressure monitor is not suitable for people suffering from serious arrhythmia

- In the event of error, an error message is displayed on the screen Er_.
- Error messages may be displayed when:
 - the systolic or diastolic pressure has not been measured (ER1 or ER2 is displayed on the screen);
 - systolic or diastolic pressure is outside the measuring range ("Hi" or "Lo" is displayed on the screen);
 - the cuff is too tight or too loose (ER3 or ER4 is displayed on the screen);
 - blood pressure is higher than 300 mmHg (ER5 is displayed on the screen);
 - inflation takes more than 180 seconds (ER6 is

displayed on the screen);

- there is an error on the system or the device (ER7 or ER8 or ERA is displayed on the screen);
- the batteries are nearly flat (18).

Technical alarm – Description

If the (systolic or diastolic) blood pressure taken is outside the range indicated in the “Technical Features” paragraph, the technical alarm displays the message “Hi” or “Lo” on the screen. In this case, consult a doctor or check that you are using the device properly. The technical alarm’s limit values are values set by the factory and can neither be changed nor deactivated. Under the IEC 60601-1-8, these limit values of the alarm have a lower priority. The technical alarm is not a locked alarm and does not need to be reset. The signal displayed on the screen switches off automatically after approximately 8 seconds.

MAINTENANCE AND STORAGE

- Clean the device and the cuff carefully, using only a slightly damp cloth.
- Do not use detergent or solvent.
- Never put either the device or the cuff in water as this might seep into the device or cuff and damage them.
- Do not put heavy objects on the device or on the cuff when they are being stored. Remove the batteries.
- No component can be maintained by user in the monitor. The circuit diagrams, component part lists, descriptions, calibration instructions, or other information which will assist the user’s appropriately qualified technical personnel to repair those parts of equipment which are designated repairably can be supplied.
- It is recommended the performance should be checked every 2 years or after repair. Please contact the service center.
- The monitor can maintain the safety and performance characteristics for a minimum of 10,000 measurements or three years.
- The cuff integrity is maintained after 1,000

open–close cycles of the closure.

- It is recommended the cuff should be disinfected 2 times every week if needed (For example, in hospital or in clinique). Wipe the inner side (the side contacts skin) of the cuff by a soft cloth squeezed after moistened with Ethyl alcohol (75-90%), then dry the cuff by airing.
- Clean the cuff after the usage of every 200 times is recommended.

TECHNICAL SPECIFICATIONS OF THE PRODUCT

Model number	LA090206
Measuring method	Measuring blood pressure on the arm, oscillometric and non-invasive
Measuring range	Pressure in the cuff 0 – 300 mmHg Systolic pressure 60 – 260 mmHg Diastolic pressure 40 – 199 mmHg Pulse 40 – 180 beats/minute
Accuracy of the indicator	systolic ± 3 mmHg diastolic ± 3 mmHg less than 60: ± 3 bpm, more than 60(incl.): $\pm 5\%$
Memory	4 x 30 memory locations
Dimensions	+/- 107 x 80 x 52mm
Weight	+/- 287g (without any batteries)
Cuff size	22 - 30cm
Admin. operating conditions	+ 10°C a + 40°C, $\leq 85\%$ relative air humidity (non-condensing)
Permissible storage and transport conditions	-20°C a +50°C, $\leq 85\%$ relative air humidity, 800 – 1 050 hPa of ambient pressure
Power supply	4 x 1.5 V AAA batteries
Battery lifespan	Approximately 200 measurements, depending on the increased blood pressure and inflation pressure
Accessories	Instructions for use, storage box
Ranking	Internal power supply, IP20, no AP or APG, continuous use, BF device
Pressure measuring system components	Pump, Valve, LCD, Cuff, Sensor

Changes could be made to the technical features without any prior warning for updating purposes.

- The Electronic Sphygmomanometer corresponds to the below standards: IEC 60601-1: Edition 3.1 2012-08/EN 60601-1:2006/A1:2013

(Medical electrical equipment -- Part 1: General requirements for basic safety and essential performance), IEC60601-1-2:2014/EN 60601-1-2:2015(Medical electrical equipment -- Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests), IEC80601-2-30:2009+AMD1:2013/ EN 80601-2-30:2010/A1:2015(Medical electrical equipment -- Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers) EN 1060-1: 1995 + A2:2009 (Non-invasive sphygmomanometers - Part 1: General requirements), EN 1060-3: 1997 + A2:2009 (Non-invasive sphygmomanometers - Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems);ISO81060-2: 2013(Non-Invasive Sphygmomanometers - Part 2: Clinical Validation Of Automated Measurement Type).

The accuracy of this blood pressure monitor has been tested correctly and its durability has been designed for long-term use. Technical measurement checks must be carried out using the appropriate methods when the device is used for medical purposes. For accurate information on how to check the accuracy of the device, you can send a request by post to the after-sales service.

ADVICE REGARDING THE DISPOSAL OF WASTE

All the packaging is composed of materials that pose no hazard for the environment and which can be disposed of at your local sorting centre to be used as secondary raw materials. The cardboard may be disposed of in a paper recycling bin. The packaging film must be taken to your local sorting and recycling centre.

When you have finished using the device, please dispose of it in an environmentally friendly way and in accordance with the law. Prior to disposal, please remove the battery and dispose of it at a collection point so that it can be recycled.

▲ Under no circumstances should spent batteries be disposed of with household waste.

EC 0197 DECLARATION DE CONFORMITY

Medical device
Category IIa.

We, LANAFORM,
hereby declare that the medical device
"ABPM-100" Lot001 LA090206 complies with European
Directive 93/42/EC concerning medical devices and
its most recent amendment 2007/47/EC.

LIMITED WARRANTY

LANAFORM guarantees this product against any material or manufacturing defect for a period of two years from the date of purchase, except in the circumstances described below.

The LANAFORM guarantee does not cover damage caused as a result of normal wear to this product. In addition, the guarantee covering this LANAFORM product does not cover damage caused by abusive or inappropriate or incorrect use, accidents, the use of unauthorized accessories, changes made to the product or any other circumstance, of whatever sort, that is outside LANAFORM's control.

LANAFORM may not be held liable for any type of circumstantial, indirect or specific damage.

All implicit guarantees relating to the suitability of the product are limited to a period of two years from the initial date of purchase as long as a copy of the proof of purchase can be supplied.

On receipt, LANAFORM will repair or replace your appliance at its discretion and will return it to you. The guarantee is only valid through LANAFORM's Service Centre. Any attempt to maintain this product by a person other than LANAFORM's Service Centre will render this guarantee void.

ELECTROMAGNETIC COMPATIBILITY INFORMATION

Table 1 - Emission

Phenomenon	Compliance	Electromagnetic environment
RF emissions	CISPR 11 Group 1, Class B	Home healthcare environment
Harmonic distortion	IEC 61000-3-2 Class A	Home healthcare environment
Voltage fluctuations and flicker	IEC 61000-3-3 Compliance	Home healthcare environment

Table 2 - Enclosure Port

Phenomenon	Basic EMC standard	Immunity test levels
		Home Healthcare Environment
Electrostatic Discharge	IEC 61000-4-2	±8 kV contact ±2kV, ±4kV, ±8kV, ±15kV air
Radiated RF EM field	IEC 61000-4-3	10V/m 80MHz-2.7GHz 80% AM at 1kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	Refer to table 3
Rated power frequency magnetic fields	IEC 61000-4-8	30A/m 50Hz or 60Hz

Table 3 – Proximity fields from RF wireless communications equipment

Test frequency (MHz)	Band (MHz)	Immunity test levels
385	380-390	Pulse modulation 18Hz, 27V/m
450	430-470	FM, ±5kHz deviation, 1kHz sine, 28V/m
710	704-787	Pulse modulation 217Hz, 9V/m
745		
780		
810		
870	800-960	Pulse modulation 18Hz, 28V/m
930		
1720		
1845	1700-1990	Pulse modulation 217Hz, 28V/m
1970		
2450		
5240	2400-2570	Pulse modulation 217Hz, 28V/m
5500	5100-5800	Pulse modulation 217Hz, 9V/m
5785		

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**ANDON HEALTH CO., LTD.**
**IHEALTHLABS EUROPE SAS**