LANAFORM

HEALTH

CN-100P

COMPRESSOR NEBULIZER

LOWER & UPPER
AIRWAYS

ADULT
CHILD

BATTERY

USB
POWERED

MMAD
3 μ m

NEBULISATION
RATE

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**DRUK-
VERSTUIVER**

**KOMPRESSI-
ONVERNEBLER**

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
EN Instruction manual

C The photographs and other images of the product in this manual and on the packaging are as accurate as possible, but do not claim to be a perfect reproduction of the product.

Thank you for choosing the LANAFORM CN-100P nebulisation kit. Combining technical sophistication and ease of use, the CN-100P nebulisation kit helps prevent diseases of the upper and lower respiratory tracts, limit their side effects and speed up recovery. The compressor propels air into the nebulisation kit. When the air penetrates the nebulisation kit, the prescribed drug is converted into an aerosol of microscopic droplets than can be easily inhaled. This product is a medical device and designed to be inhaled at home. Drugs must only be inhaled after taking medical advice.

OPERATION OF THE NEBULISATION KIT

The drug that is pumped through the drug channel is mixed with compressed air generated by a compressor pump. The compressed air mixed with the drug is converted into fine particles and vapourised as it collides with the deflector.

SYMBOLS USED

The following symbols are used in the instructions for use.

	Warning	This symbol alerts you to potential hazards that can injure you or damage your health.
	N.B.	This symbol alerts you to potential damage that may be caused to the device or to an accessory.
	Note	This symbol indicates important information.

TITLE AND DESCRIPTION OF DISPLAY SYMBOLS ON THE SHEET DESCRIBING THE PRODUCT'S FEATURES

The following symbols are used on the packaging and on the data plate of the device and its accessories.

Identification	Definition of the symbol
	This EC symbol demonstrates that the device complies with the basic requirements of directive 93/42/EEC on medical devices
	European representative: Shanghai international Holding Corp. GmbH (Europe) Address: Eiffeustrasse 80, 20537 Hamburg Germany
	Manufacturer's details: Shenzhen Kingyield Technology Co., Ltd. A5 Bldg, FuZhong Industrial Zone, Xia Shi Wei road, FuLong district, BaoAn Shenzhen, China
	Batch number
	BF device
	Instructions for use. Follow the instructions on how to use the device
	Protective rating for liquids and solids
	The device must be disposed of in accordance with the European WEEE (Waste Electrical and Electronic equipment) Directive on electrical and electronic devices.

PLEASE READ ALL THE INSTRUCTIONS BEFORE USING YOUR DEVICE, ESPECIALLY THIS BASIC SAFETY INFORMATION

Read all the information in the instructions for use and in the documentation included in the packaging carefully before using the device.

When an electrical device is in use, and particularly when children are present, basic safety precautions must always be taken, specifically the following:

- When using this product, please make sure you comply with all the items listed below.
- Follow your doctor's instructions and those of health-care staff as regards the nature of the drug, its prescription instructions and dosage.
- If you experience any unusual sensations when using the device, stop using it immediately and consult your doctor.
- Do not use water on its own in the nebuliser for inhaling.
- Clean and disinfect the nebulisation kit, the child's and adult's mask and the mouth piece before using it for the first time if the device has not been used for a lengthy period or if several people use the same device.
- Make sure you wash the components after use and ensure they are thoroughly disinfected, dried and stored in a clean place.
- Keep the device and child's accessory out of the reach of babies and children not being supervised. The device is likely to contain small parts that can be swallowed.
- Do not put the air tube away if there is any condensation or any of the drug left inside it.
- Do not leave cleaning solution in the nebuliser's components. Rinse them in hot tap water after disinfecting them.
- Do not use or keep the device in places where it might be exposed to harmful vapours or volatile substances.
- Do not use or keep the device in places where it might be exposed to flammable gases.
- Do not cover the compressor with any cover, towel or anything else when it is in use.
- Always dispose of any remaining drug after use and use a new dose of the drug each time.

Risk of electric shock

- Do not use the compressor or USB cable if they are wet and do not connect them to a power outlet or other devices with wet hands.
- The compressor is not water-proof. Do not spill water or any other liquids on these components. If a liquid is spilled onto these components, disconnect the USB cable immediately and soak up the liquid with a gauze pad or other absorbent fabric.
- Do not immerse the compressor in water or any other liquid.

- Do not use the device or store it in damp places such as a bathroom.
- Do not use the device if the USB cable is damaged.
- Keep the USB cable away from hot surfaces.
- Read and note the "Important information about electromagnetic compatibility (EMC)" displayed in the Technical Data section.
- This device must not be left unsupervised if used by, or close to, children or the infirm.
- Check that the components have been fixed correctly.
- Make sure that the spray cap has been installed correctly prior to use.
- Check that the air filter has been fixed correctly.
- Check that the air filter is clean. If the air filter has changed colour or been in regular use for over 60 days, replace it with a brand new air filter.
- Do not spill any liquid or any of the drug on the compressor.
- Do not tilt the nebulisation kit at an angle of more than 45° in any direction and do not shake it while it is running.
- Do not use or store the device if there are folds in the air tube.
- Use only original nebulisation components, air tube, air filter and air filter cover.
- Do not add more than 10ml of the drug to the drug reservoir.
- Do not carry or remove the nebulisation kit when the drug reservoir contains a drug.
- Do not leave the device unsupervised in the presence of young children or individuals who cannot give their consent.
- Do not subject the device or its components to violent shocks such as, for example, dropping the device onto the floor.
- Do not twist the spray cap or push a pin or pointed object through the cannula of the drug reservoir.
- Do not insert your fingers or any objects in the compressor.
- Do not dismantle or attempt to repair the compressor or USB cable.
- Do not leave the device or its components in places where the temperature is beyond the range of those indicated in the technical data of these instructions for use.
- Do not leave the device or its components in places where they will be exposed to extreme temperatures or changes in humidity, such as in a vehicle at the height of summer or in places where they will be exposed to direct sunlight.
- Do not block the air filter cover.
- When you disinfect components by boiling them, make sure the container is not dry-heated.
- Never use the device when you are going to sleep or are drowsy.
- Do not obstruct the slot between the cartridge and the air inlet of the inhalation system.
- Do not use the microwave oven, dishwasher drier or hairdryer to dry the device or its components.
- Authorised for human use only.
- Continuous use of the device may reduce its lifespan. Do not use this device for more than 15 minutes at a time and wait for at least 40 minutes before using it again.
- The main unit may get hot when the device is in use.
- Do not touch the main unit except when necessary, for example to switch the device off during nebulisation.
- To prevent any of the drug remaining on your face, make sure you cleanse your face after taking off the mask.

Risk of electric shock


- Always remove the batteries or always disconnect the USB cable after use and before cleaning.
- Any changes or modifications that have not been approved by LANAFORM will cancel the user's warranty.

GENERAL SAFETY PRECAUTIONS

- Examine the device and the components before each use and make sure that no anomalies are found. Check the following points in particular:
 - The spray cap, cannula, child's/adult's mask and the air tube must not be damaged.
 - The cannula must not be obstructed.
 - The compressor must run normally.
- The compressor pump emits sounds and vibrations when this device is in use. Compressed air emission from the nebulisation kit also generates noise. This noise is normal and is not the sign of a malfunction.
- Use the device only for the intended purposes. Do not use it for any other purpose.
- Do not use the device at temperatures above +40°C and below 5°C.
- Check that the air tube is correctly attached to the compressor and nebulisation kit and that it does not become loose. When you insert the air tube plug into the air tube connectors, screw it loosely to prevent the air tube becoming disconnected when in use.
- Remove the batteries if the device is not used for more than 3 months.
- Read and take note of the "How to dispose of this product" procedure in the Technical Data section when you discard the device and the old optional accessories or components.
- Look after these instructions carefully so that you can refer to them at a later date.

**DESCRIPTION OF THE DEVICE
AND ITS ACCESSORIES**

- 1 Compressor
- 2 Nebulisation kit bracket
- 3 Compressed air tube connector
- 4 Operating light indicator
- 5 Stop/start button
- 6 USB connector
- 7 Air filter cover
- 8 Battery cover
- 9 Compressed air tube
- 10 Nebulisation kit
- 11 Cartridge
- 12 -
- 12A Top section of the inhalation system
- 12B Air inlet of the inhalation system
- 13 Spray cap
- 14 Drug reservoir
- 14A Cannula
- 14B Air tube connector
- 15 Children's mask
- 16 Adult's mask
- 17 Mouth piece
- 18 Replacement filter
- 19 Storage bag

PREPARING THE NEBULISER FOR USE

⚠ Clean and disinfect the nebulisation kit, the mask and mouth piece before using it for the first time after buying it. Clean and disinfect the nebulisation kit, the mask and the mouth piece before using it if the device has not been used for a lengthy period or if several people use the same device.


For cleaning and disinfecting instructions, consult the appropriate chapter of these instructions for use.

ⓘ Make sure that the air filter is changed regularly. Consult the section on "Changing the air filter".

1. Remove the nebulisation kit cartridge.

2. Remove the top section of the drug reservoir's inhalation system.

- a) Turn the top part of the inhalation system in an anti-clockwise direction.
- b) Lift the top section of the inhalation system to remove it from the drug reservoir.


3. Add the appropriate volume of the prescribed drug to the drug reservoir.

⚠ The capacity of the drug reservoir is between 2 and 10ml.

ⓘ The scale markings that appear on the nebulisation kit are shown for reference only. Please use the scale markings on your syringe or bottle for accurately measuring the drug.


4. Check that the spray cap is correctly fixed to the top part of the inhalation system.


5. Fix the top section of the inhalation system to the drug reservoir.


- a) Place the top section on the drug reservoir, as illustrated.
- b) Turn the top part of the inhalation system in a clockwise direction and align it with the drug reservoir until they click into place, as illustrated.

6. Fix the cartridge to the top section of the inhalation system.


7. Place the batteries (4 AA 1.5V – not supplied) in the container provided for this purpose or connect the USB cable supplied with the nebulisation kit and an adapter of 1.0A minimum (not supplied).


8. Fix the accessory required to the top section of the inhalation system.

- a) Fix the child's or adult's mask to the top section of the inhalation system.
- b) Fix the mouth piece to the top section of the inhalation system.


9. Fix the air tube.

- a) Loosely screw the air tube bracket and insert it firmly into the air tube connector on the compressor.
- b) Loosely screw the air tube bracket and insert it firmly into the air tube connector on the drug reservoir.

ⓘ Check that the air tube is correctly attached to the compressor and nebulisation kit and that it does not become loose. When you insert the air tube plug into the air tube connectors, screw it loosely to prevent the air tube becoming disconnected when in use.

Make sure that you do not spill the drug when you are attaching the air tube. Keep the nebulisation kit upright. Use the nebulisation kit bracket as a temporary support for the nebulisation kit.


10. Description of the compressor's status

Once you press the stop/start button, the compressor starts running and the operating light indicator lights up. The compressor will run for 15 minutes.

If you press the stop/start button a second time, the compressor stops and the light indicator goes out.

If the batteries are flat, the compressor stops and the light indicator flashes for 2 minutes, indicating that the batteries need to be changed.

If possible, choose to run the device using the USB cable rather than batteries.

Indicator Light

ADMINISTERING THE DRUG

1. Hold the nebulisation kit as shown on the right.

⚠ Do not tilt the nebulisation kit at an angle of more than 45° in any direction. The drug may run into your mouth or nebulisation might not be effective.

2. Turn the mains switch to the ON position (I). The compressor switches on and nebulisation starts.

To stop nebulisation Turn the mains switch to the OFF position (II).


⚠ Do not cover the compressor with any cover or towel etc. when it is in use.

Do not block the air filter cover.

3. Inhale the drug as instructed by your doctor and/or your therapist.

a) Using the adult's/child's mask


Attach the mask so that it covers your nose and mouth and inhale the drug. Breathe out through the mask.

b) Using the mouth piece

Insert the mouth piece into your mouth and inhale the drug by breathing normally.

Breathe out normally through the mouth piece.



4. Switch off the device at the end of treatment. Check that there is no condensation or moisture in the air tube and remove the batteries or disconnect the USB cable and the IA minimum adapter (not supplied) that has been used.

⚠ Some condensation may form in the air tube. Do not put the air tube away if there is any condensation or moisture left inside it. It might be infected by bacteria.

If there is moisture or liquid still left in the air tube, remove the moisture as follows.

- Check that the air tube is still connected to the air intake duct on the compressor.
- Check that the compressor is connected to the appropriate USB cable or that the batteries are inserted correctly in their compartment.
- Disconnect the air tube from the nebulisation kit.
- Start the compressor and pump the air through the air tube to expel the moisture.

DAILY CLEANING AND DISINFECTING ROUTINE

A. CLEANING


If you follow the cleaning instructions after each use, you will prevent the drug left behind in the reservoir from drying, thus allowing the device to nebulise efficiently as well as preventing any contamination.

⚠ Rinse the nebulisation kit, children's accessory, mask or mouth piece thoroughly with hot tap water each time it is used and disinfect it after the last treatment of the day.


- Check that the mains switch is in the OFF position ().
- Disconnect the USB cable from the compressor's connector and/or take the batteries out of their compartment.
- Disconnect the nebulisation kit from the air tube and dismantle it.

- Dispose of any of the drug left behind in the drug reservoir.
- Wash all the nebulisation components with hot water and mild detergent and rinse with hot tap water.
- Having cleaned and disinfected the components, rinse them thoroughly in clean water and allow to dry in the air in a clean place.

- ① Do not start the device up until the clean components are completely dry.
- 7 If the compressor or outer section of the air tube is dirty, clean them with a soft cloth moistened with water or mild detergent.

- ⚠ The compressor is not water-proof.
- 8 Reassemble the nebulisation kit. Put it in the storage bag.

B. DISINFECTING THE NEBULISATION COMPONENTS

Always disinfect the drug reservoir, mouth piece and other nebulising components after the last treatment of the day. To disinfect the components, follow one of the methods described below:

- B.1. Use a commercially available disinfectant. Follow the manufacturer's instructions.
 - 1 Immerse the components in the cleaning solution for the period indicated.
 - 2 Remove the components and discard the solution.
 - 3 Rinse the components thoroughly in hot tap water, remove the surplus and allow to dry in the air in a clean place.
- B.2. You can boil the components for 15 to 20 minutes. Carefully remove the components after boiling, remove excess water and allow to dry in the air in a clean place.

C. HANDLING THE SPRAY CAP

The spray cap is an important component for spraying the drug. Please take the following precautions when handling it.

- ⚠ Wash the spray cap after each use.
- Do not use the brush, pin, etc. to clean the components.
- If you disinfect the components in boiling water, make sure you boil them in large volumes of water.
- Do not boil the spray cap with items other than the accessories corresponding to the nebulisation system.
- Do not use the microwave oven, dishwasher drier or hairdryer to dry the device or its components.

MAINTENANCE AND STORAGE

- ⚠ To keep your device in the best possible condition and in order not to damage it, proceed as follows:
 - Do not leave the device unsupervised in the presence of young children or individuals who cannot give their consent.
 - Do not subject the device or its components to violent shocks such as, for example, dropping the device onto the ground.

- Do not carry or withdraw the nebulisation kit when the drug reservoir contains a drug.
 - Do not store the device in a place that is subject to extreme hot or cold temperatures, high humidity or in direct sunlight.
 - Do not bend or fold the air tube.
 - Dispose of the device and its optional used accessories or components according to the local laws and regulations.
 - Never use benzene, thinner or flammable chemical product to clean the device.
- Always carry the device in its storage bag and store it in the bag if you do not anticipate using it for a long time.

CHANGING THE AIR FILTER

If the air filter has changed colour or been in regular use for over 60 days, replace it with a brand new air filter.

- 1 Remove the compressor's air filter cover.
- 2 Change the air filter.
- 3 Put the air filter cover back in situ.

- ①
 - Use only air filters supplied with the device, which are designed specifically for this device. Do not operate the unit without a filter.
 - The air filters can be inserted in any direction.
 - Check that the air filters are clean and free of dust before inserting them.
 - Do not attempt to wash or dry the filter. If the air filter is wet, replace it. Damp air filters may cause obstructions.
 - To prevent the air filter cover becoming clogged, wash it regularly. Do not let the air filter cover boil. Remove the air filter cover, wash it and make sure it is completely dry and attached correctly.
 - Please contact your LANAFORM reseller to obtain additional air filters or if the components supplied with the device are damaged or missing.

TROUBLE-SHOOTING

Check the following points if your unit malfunctions when running. You can also consult the pages of these instructions for use for detailed instructions.

Symptoms	Cause	Solution
Nothing happens when the mains switch is in the ON position.	Have the batteries been inserted into their compartment or has the USB cable been connected to the compressor's connector?	Check that the plug has been inserted correctly into the electric socket.
Nebulisation stops automatically after 10 seconds when the device is on.	The batteries are flat.	Change the batteries or connect the USB cable after removing the batteries from their compartment.

No nebulisation or nebulisation speed slow when the device is on.	Does the drug reservoir contain any of the drug?	Add the appropriate volume of the drug to the drug reservoir.
	Does the drug reservoir contain too much or not enough of the drug?	
	Is the spray cap missing or not assembled correctly?	Attach the spray cap correctly.
	Has the nebulisation kit been correctly assembled?	Assemble the nebulisation kit correctly.
	Is the cannula blocked?	Check that the cannula is not obstructed.
	Is the nebulisation kit tilted at too steep an angle?	Make sure that the nebulisation kit is not tilting at an angle of more than 45 degrees.
	Is the air tube connected correctly?	Check that the air tube is correctly attached to the compressor and nebulisation kit.
	Is the air tube folded or damaged?	Check that the air tube is not twisted.
	Is the air tube blocked?	Check that there air tube is clear of obstructions.
	Is the air filter dirty?	Replace the air filter with a brand new one.
The device is unusually noisy.	Has the air filter cover been correctly attached?	Attach the air filter cover correctly?
The device is very hot.	Is the compressor covered?	Do not cover the compressor with anything when in use.

ⓘ

If the solution suggested does not resolve the problem, do not try to repair the device since no component of the unit can be repaired by the user.

Return the unit to a retail outlet or approved LANAFORM distributor.

TECHNICAL INFORMATION

- Name: LANAFORM
- Model: CN-100P LA120420
- Type: Compression-driven nebuliser
- Nominal value (compressor with nebuliser): 1.5V x 4 AA batteries or USB 5 V 1 A
- Operating conditions: Operates intermittently: 20 mins. on/40 mins. off
- Temperature/Moisture content when running: +5°C to +40°C/15% at 85% HR
- Temperature/Moisture content when stored: -25°C to +35°C with relative humidity below 90°C with no condensation and between 35°C and 70°C with water vapour up to 50hPa
- Atmospheric pressure: HR 700 - 1,060 hPa
- Weight: Approximately 180g (compressor only)
- Dimensions: Approximately 115 (L) x 77 (H) x 42 (P) mm
- Contents: Compressor, nebulisation kit, tube à air (PVC, 100cm), replacement air filters (5 units), mouth piece, adult's mask (PVC), children's mask (PVC), USB cable, storage bag, instructions for use.
- Operating time: 15 minutes on/40 minutes off.
- Built-in timer: The compressor runs for 15 minutes once the on/off button has been pressed and then stops.
- Compression output: 2-10/min
- Compression pressure: 360-550mmHg

ⓘ

- May be changed without prior notice.
- The device may not run if the temperature and voltage are different from those defined in the technical features.
- Do not use the device in places where it might be exposed to flammable gases.
- This device complies with the provisions of European Community Directive 93/42/EEC (Directive on medical devices) and European standard EN13544-1:2007, Oxygen Therapy Equipment - Part 1: Nebulisation systems and their components.

TECHNICAL INFORMATION ON THE COMPRESSOR WITH NEBULISATION KIT

Volume of the appropriate drug:	2ml minimum - 10ml maximum
Drug reservoir capacity:	10ml maximum
Size of the particles:	**MMAD Approximately 3.9 µm MMAD = mass median aerodynamic diameter
Nebulisation rate:	Approximately 0.3ml/min (depending on weight loss)
Acoustics:	Approximately 58dB at a distance of 1 metre
Aerosol outlet:	0,3627 ml (1 % M/V) NaF
Aerosol output:	0,045 ml/min (1 % M/V) NaF/min

ⓘ

Performance may vary with certain drugs such as suspensions or drugs with high viscosity. Consult the drug supplier's data sheet for further details.

** Measurements taken independently according to standard EN 13544-1.

IMPORTANT INFORMATION ABOUT ELECTROMAGNETIC COMPATIBILITY (EMC)

With the increase in the number of electronic devices such as the PC and mobile (cell) phones, medical devices can be subject when in use to electromagnetic interference released by other devices. Electromagnetic interference can disturb the operation of the medical device and create a potentially dangerous situation.

Medical devices must not interfere with other devices either.

In order to regulate the requirements relating to the EMC (electromagnetic compatibility) with the aim of preventing any dangerous situation arising caused by the product, the EN 60601-1-2:2007 standard has been implemented. This standard defines the degrees of immunity to electromagnetic interference as well as the maximum levels of electromagnetic emissions for medical devices.

This medical device complies with this EN 60601-1-2:2007 standard both in terms of immunity and emissions.

It is nevertheless important to take special precautions:

Do not use mobile (cell) phones and other devices that generate powerful electrical or electromagnetic fields close to the device. This might disturb the operation of the medical device and create a potentially dangerous situation. Keeping a minimum distance of 7m away is recommended. Check that the device is running properly if it is closer than this.

♻️ ADVICE REGARDING THE DISPOSAL OF WASTE

All the packaging is composed of materials that pose no hazard for the environment and which can be disposed of at your local sorting centre to be used as secondary raw materials. The cardboard may be disposed of in a paper recycling bin. The packaging film must be taken to your local sorting and recycling centre.

When you have finished using the device, please dispose of it in an environmentally friendly way and in accordance with the law. Prior to disposal, please remove the battery and dispose of it at a collection point so that it can be recycled.

Please note: under no circumstances should spent batteries be disposed of with household waste.

EC 0197 DECLARATION OF CONFORMITY

Medical device
Category IIa.
We, LANAFORM,

hereby declare that the medical device « CN-100P » Lot001 LA0900420 complies with European Directive 93/42/EC concerning medical devices and its most recent amendment 2007/47/EC.

LIMITED WARRANTY

LANAFORM guarantees this product against any material or manufacturing defect for a period of two years from the date of purchase, except in the circumstances described below.

The LANAFORM guarantee does not cover damage caused as a result of normal wear to this product. In addition, the guarantee covering this LANAFORM product does not cover damage caused by abusive or inappropriate or incorrect use, accidents, the use of unauthorized accessories, changes made to the product or any other circumstance, of whatever sort, that is outside LANAFORM's control.

LANAFORM may not be held liable for any type of circumstantial, indirect or specific damage.

All implicit guarantees relating to the suitability of the product are limited to a period of two years from the initial date of purchase as long as a copy of the proof of purchase can be supplied.

On receipt, LANAFORM will repair or replace your appliance at its discretion and will return it to you. The guarantee is only valid through LANAFORM's Service Centre. Any attempt to maintain this product by a person other than LANAFORM's Service Centre will render this guarantee void.

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