

Cell Roller







What does "squeeze-roll" mean?

The squeeze-roll technique is an ancestral method that has proven its worth. In concrete terms, this technique involves pinching the skin between the fingers and pulling it upwards with the "fat layer".

The main purpose of squeeze-roll is to **break** up the fatty deposits in the skin, to help get rid of the orange peel effect.

As well as its anti-cellulite properties, the miraculous squeeze-roll technique also improves blood circulation and lymph flow, to deal with the symptoms of heavy and puffy legs.

Remember, this technique goes hand in hand with a healthy lifestyle.





In practice

The massage should be performed daily to see rapid results. Use the Cell Roller for **5 minutes per area.**The device can be used with a massage cream to make the session a more enjoyable experience. Most of all, **the skin will be nourished and hydrated.**

The Cell Roller is **easy to use** and its **limited size** is an asset where there is a lack of space and storage. **It works** well on the thighs, buttocks, belly, waist and arms.





Cell Roller THE ANTI-CELLULITE MASSAGE DEVICE

LA110223

Features

- Anti-cellulite and anti-orange peel effect
- Squeeze-roll massage for significantly firmer and smoother skin
- Encourages blood flow and stimulates lymphatic circulation
- Suitable for the thighs, arms, belly and buttocks
- High-performance motor for a powerful squeeze-roll massage.
- 2 levels of massage intensity
- 3x2 rotating massage rollers
- Ergonomically designed for easy and practical use
- Easy to clean, courtesy of the detachable massage rollers







Outer packaging dimensions

 $41.0 \times 30.36 \times 29.6 \text{ cm}$



Outer packaging weight

6.5 kg



Colour Boxes/ Outer packaging

12



Colour Box Dimensions

 $13.8 \times 19.6 \times 9.6 \text{ cm}$



Colour Box Weight

0.5 kg



Voltage

12 V



Home Health Experience

Rue de la Légende 55 4141 Louveigné, Belgium

Tel.: +32(0)43609291 Fax: +32(0)43609723

www.lanaform.com info@lanaform.com







