

LANAFORM

Hydrojet

**EFFECTIVELY REMOVES
DENTAL PLAQUE**

AND MAKES IT EASIER TO CLEAN DENTAL APPLIANCES



ORAL IRRIGATOR

Improve your gum health and get a gentle and targeted cleaning of your teeth

3 specific modes for appropriate treatment

With its **3 cleaning modes**, the Hydrojet not only cleans the spaces between your teeth by gently dislodging food residues, but also by removing dental plaque and bacteria.

The oral irrigator also means that you can massage your gums gently while relieving any bleeding.

Soft




Suitable for gently massaging the gums, it gently stimulates problem areas by aiming the cannula at the affected area.

Normal

Intended for treatments that focus on periodontal pockets, it removes bacteria and dental plaque.

Pulse

This mode is intended for treatments that focus on the teeth to remove bacteria as well as tartar.





▲ **Cleans your whole mouth in just 1 minute**

With a 230ml reservoir, the hydrojet cleans your whole mouth in just 1 minute! Position the guide tip, move the cannula and repeat the operation to clean your mouth completely.



Wireless technology and easy to use

With its rechargeable battery, the hydrojet has a battery life of 45 minutes. Its ergonomic handle is designed to be easy to manoeuvre in small spaces.




The guide tip makes sure the positioning is correct

Guide tip and fine cannula that can be tilted 360° make positioning easy. Simply slide the tip over the teeth to clean the space between each tooth or gum.




Treatment for each member of the family

The 4 360° rotating cannulas mean that it can be used by each member of the family.



Features

Modes	3
Rechargeable battery	● (1800 mAh)
Battery life	45 min.
Water tightness index	IPX7
Contents	230 ml
Accessories	4 360° rotating cannulas included
Speed	1,200 pulses/min.



Master Carton Dimensions
52.5 × 40.0 × 38.7 cm



Master Carton Weight
15.75 kg



Colour Boxes / Master Carton
30



Colour Box Dimensions
12.0 × 7.5 × 25.7 cm



Colour Box Weight
0.480 kg



Voltage
DC 5V === 1A

LANAFORM

Home Health Experience

Rue de la Légende 55
4141 Louveigné, Belgium

Tel. : +32 (0)4 360 92 91

Fax : +32 (0)4 360 97 23

www.lanaform.com
info@lanaform.com

