

LANAFORM

Home Health Experience

Skin Mass

ANTI-CELLULITE MASSAGE DEVICE



* Squeeze-roll massage * Reduces the orange peel effect * Smoother, firmer skin





Le Skin Mass **The solution against cellulite**

Skin Mass, the massage device that reduces cellulite and redefines the figure.


The combination of the squeeze-roll massage technique and the variable intensity suction system allows Skin Mass to reduce cellulite while firming the body contours and improving the skin's elasticity.



What is squeeze-roll massage?

For many years now, the squeeze-roll massage technique has been a cult method in the field of slimming. This technique, known to effectively "break" cellulite, promotes lipolysis, i.e. the destruction of fat cells, while reducing the "orange peel" appearance of the skin.

The result? **Smoother, firmer skin.**



Squeeze-roll: how does it work?

Squeeze-roll (or palpate and roll) massage involves pinching the skin where the cellulite is, then rolling the skin over itself.

This massage method promotes lipolysis, i.e. the destruction of fat cells that cause cellulite. Blood circulation is activated, the skin becomes more supple and toned again, and the orange peel appearance is effectively reduced.

Skin Mass, the solution against cellulite

Skin Mass is a massage device that imitates the squeeze-roll technique. Thanks to its two rollers, it reproduces the consecutive double pinching of the skin in order to decompartmentalise cellulite in a more intense and rapid way than a manual massage.

Its variable intensity suction function allows it to be adapted to the area to be treated and to each type of skin, whether it is sensitive or not.


Skin Mass has three modes that you can perform one after the other as a routine or separately.



3 MODES

Roller rotation


Prepare your skin



1

Roller rotation + Suction


Breaks down cellulite



2

Roller rotation + Suction + Vibration

Smoothes and drains the skin



3


In practice

In order to guarantee fast visible results, it is recommended to use the Skin Mass daily on the areas to be treated.

Skin Mass is easy to use and allows you to perform a complete and professional treatment of your skin from your bathroom.

It can be used* on all areas that are prone to cellulite:

- Thighs (inner and outer)
- Buttocks
- Legs
- Knees
- Belly
- Arms



The user guide for the Skin Mass device is presented in a clean, white layout. At the top, the brand name 'LANAFORM' is displayed in a blue box. Below this, the title 'How to use Skin Mass' is followed by a list of languages in which the device is available: English, Français, Nederlands, Deutsch, Español, Italiano, Slovenčina, Slovenščina, Hrvatski, Русский Язык, and Български. A QR code is provided for more information, with the URL www.bit.ly/skinmass listed below it. A warning section with a triangle icon states: 'Warning: If you experience any discomfort, stop using the Skin Mass immediately and consult your doctor.' The guide is divided into three main sections, each with a small diagram showing the device's placement and movement on the body: 'Use on the belly and waist', 'Use on the legs', and 'Use on the arms'. The 'Use on the arms' section includes the instruction: 'Perform slow up-down movements of your arms.' A small image of the device is shown at the bottom right of the guide.



* A user guide with all the best practices is provided with the device.



Skin Mass

ANTI-CELLULITE MASSAGE DEVICE

LA11022002

Features

- Anti-cellulite and anti-orange peel effect
- Squeeze-roll massage for firmer, smoother skin
- Encourages blood flow and stimulates lymphatic circulation
- Suitable for the thighs, arms, belly and buttocks
- High-performance motor for a powerful squeeze-roll massage.
- Variable massage intensity
- 2 massage rollers
- 3 massage modes
- Elastic and adjustable handle
- Ergonomically designed for easy and practical use





Outer packaging dimensions
49.0 × 24.5 × 50.5 cm



Outer packaging weight
10.10 kg



Colour Boxes/ Outer packaging
12



Colour Box Dimensions
15.3 × 21.5 × 12 cm



Colour Box Weight
0.8 kg

LANAFORM

Home Health Experience

Rue de la Légende 55
4141 Louveigné, Belgique

Tel. : +32(0)4 360 92 91

Fax: +32(0)4 360 97 23

www.lanaform.com

info@lanaform.com

